


KaVo Turbines, straight and contra-angle handpieces.

Start safe with
affordable
KaVo quality.

KaVo **smart** series



KaVo. Dental Excellence.

KaVo Turbines, straight and contra-angle handpieces.

Durable, reliable quality made by KaVo.



Your KaVo quality benefit:

One hundred years of experience and more than **five million KaVo handpieces** sold reflect Dental Excellence throughout the world. Just four of many proofs of KaVo quality:

- > In our production, we place the highest emphasis on precision and have the strictest in-house quality standards, including: **clean room assembly in Biberach.**
- > While up to 70% of the handpieces made by other manufacturers require a major repair within the first four years, this only applies **to 30%** of KaVo handpieces. This puts KaVo on top in this market (source: GfK ZaBus survey March 2009).
- > The **repeated purchase rate of 98%** reflects the **extreme satisfaction** level of dentists with the quality of our products.
- > **Periodical certification** ensures consistently high quality. KaVo is certified according to EN ISO 9001 and EN ISO 13485.

Your KaVo system solution benefit:

Benefit from your KaVo System, with perfectly matched solutions in terms of **efficiency and protection of investment:**

- > **Avoid costly new purchases** and simply upgrade your handpiece bulb for the new KaVo MULTI LED.
- > **The long-term availability of replacement** parts significantly increases the protection of your investment.
- > With all new developments, we place much emphasis on matching system solutions allowing you to retro-fit your system at a later time **in a cost-efficient manner.**

Your KaVo innovation benefits:

Through innovation originating from passion and **more than 2,200 patents** and registered designs both domestic and abroad, we are contributing to your everyday success. Just a few successful examples:

- > **KaVo COMFORTdrive**, the first contra-angle handpiece with an integrated electric motor combines power and lightweight design.
- > The **exchangeable KaVo Spray Microfilter** optimises spray cooling and reduces service costs.
- > The **turbine's automatic pressure regulating valve** ensures consistently good performance whilst maintaining powerful torque.
- > The **patented head and knee angles of 100° and 19°**, respectively, of KaVo handpieces ensure optimal access to the preparation area.
- > **SONICflex Quick-Coupling** accelerates and simplifies the exchange of scaler tips.

Your KaVo safety benefits:

KaVo provides the foundations **for optimal treatment safety** and protects the dentist, assistant, and patient against infections during every treatment. How this is done by KaVo is evident from some practical examples:

- > The **collet retention** force of up to 30 N in every product line far exceeds any standard requirements and ensures optimal support for your burs in the long-term.
- > **Innovative hygiene solutions**, such as the built-in non-return valve in the GENTLEsilence series and the easily cleaned and sterilisable surfaces, protect all individuals involved against cross-contamination.

KaVo smart series

1. Made in Germany.

2. Premium materials.

3. Entry level price for the world of KaVo.



More than **5** million handpieces sold throughout the world



KaVo. Dental Excellence.

KaVo smart series

Entry level price for the world of KaVo – the handpieces of the KaVo Smart Series.

If you aim to reach the heights, you need solid foundations. The KaVo Smart Series offers the ideal prerequisites. Enter the world of Dental Excellence and experience KaVo quality and reliability, at an unbeatable price.

KaVo Smart Series – Start safe with affordable KaVo quality.




1. Made in Germany.

> **Proven quality**, that has made KaVo into the worldwide market leader

> **100 years of experience** goes into each product

> **High quality preparation**, such as one-piece handpiece casings




2. Solid, basic configuration.

> **High torque**, thanks to 16 Watt power output, for rapid, efficient preparation

> **KaVo safety standards**, such as a chuck retention force of up to 30 N


> **Always the right speed** for each individual clinical indication, thanks to various gearbox ratios



3. Economical entry into the world of KaVo.

> **Low price** for entry into the KaVo world of dental excellence

Turbines with MULTIflex connection




SMARTtorque S619 L/C / S615 L Mini

- > Light and non light version
- > **Head versions: standard** (h: 14.5 mm) or **mini** (h: 12.3 mm)
- > **Reliable 16 Watt power**, for efficient and safe preparation
- > **Different model / price versions**
- > **High safety**, with chuck retention force up to 30 N
- > **4-port spray**, for optimum cooling effect
- > Enhanced, **anti-reflux effect**
- > New long-life, **ceramic bearings**

For order details, consult latest price list

Mat. no.: S619 L - 1.008.1641
S619 C - 1.008.1642
S615 L - 1.008.1643

Turbines with fixed connection



SMARTtorque S609 C / S605 C Mini

- > **4-hole standard connection** (Midwest)
- > **Head versions: standard** (h: 14.5 mm) or **mini** (h: 12.3 mm)
- > **Reliable 16 Watt power**, for efficient and safe preparation
- > **Different model / price versions**
- > **High safety**, with chuck retention force up to 30 N
- > **4-port spray**, for optimum cooling effect
- > Enhanced, **anti-reflux effect**
- > New long-life, **ceramic bearings**

For order details, consult latest price list

Mat. no.: S609 C - 1.008.1644
S605 C - 1.008.1645

SMARTtorque S608 C

- > **2-hole standard connection** (Borden)
- > **Head version: standard** (h: 14.5 mm)
- > **Reliable 16 Watt power**, for efficient and safe preparation
- > **High safety**, with chuck retention force up to 30 N
- > **4-port spray**, for optimum cooling effect
- > Enhanced **anti-reflux effect**
- > **KaVo Push-button chuck system**, for rapid bur changing
- > **Ceramic ball bearing**

Mat. no. 3.001.6560

12 months warranty.

MULTIflex couplings



MULTIflex LUX coupling 460 LED

- > With **KaVo MULTI LED**
- > **More flexibility**, due to a connection that is compatible with many handpieces
- > **Increased efficiency**, as fast handpiece interchange is possible

For all KaVo turbines with MULTIflex connection, SONICflex airscaler, PROPHYflex polisher, CORONAflex, all INTRAflex LUX motor-handpieces, airmotor 181 M, RONDOflex

Mat. no. 1.007.3201




MULTIflex LUX coupling 460 LE

- > With the proven, **KaVo high-pressure, halogen bulb**
- > **More flexibility**, due to a connection that is compatible with many handpieces
- > **Increased efficiency**, as fast handpiece interchange is possible

For all KaVo turbines with MULTIflex connection, SONICflex airscaler, PROPHYflex polisher, CORONAflex, all INTRAflex LUX motor-handpieces, airmotor 181 M, RONDOflex

Mat. no. 1.001.7599



MULTIflex coupling 460 E

- > **More flexibility**, due to a connection that is compatible with many handpieces
- > **Increased efficiency**, as fast handpiece interchange is possible

For all KaVo turbines with MULTIflex connection, SONICflex airscaler, PROPHYflex polisher, CORONAflex, all INTRAflex LUX motor-handpieces, airmotor 181 M, RONDOflex

Mat. no. 1.001.7600



SMARTmatic Series – Quality made in Germany. Variety made for you.

The new SMARTmatic series unites the legendary, Made in Germany KaVo quality with smart features for optimal ease of use and clinical excellence. As a successor to the established INTRAmatic E / ES series, our SMARTmatic series fits alongside the MASTER and EXPERT series to offer you a new attractive entry into the world of KaVo.

Whichever instrument you choose from the KaVo SMARTmatic series, it will amaze you with its many benefits that other instruments just cannot offer.



100% Made in Germany

All instruments are developed and produced exclusively in Germany.



Premium materials

Rust-free stainless steel for long-lasting durability and low running costs.



Universally suitable

Usable with all dental chairs with air and electric motors (ISO 3694) and internal or external spray.



Smaller head

Compact dimensions for the best view and precise treatment.



SMARTmatic S20 S

- > Transmission ratio 1:1
 - > For **contra-angle handpiece burs with 2.35mm**
 - > **Internal spray transmission** to the Sprayclip
- Mat. no. 1.011.6752
also available without Spray:
Mat. no. 1.011.6750



SMARTmatic S10 S

- > Transmission ratio 1:1
 - > For **straight and contra-angle handpiece burs with 2.35mm**
 - > **Internal spray transmission** to the Sprayclip
- Mat. no. 1.011.6732
also available without Spray:
Mat. no. 1.011.6730



SMARTmatic ENDO S81

- > Transmission ratio 8:1
 - > For **endodontic files with shaft Ø 2.35 mm**
- Mat. no. 1.011.6780



SMARTmatic ENDO S321

- > Transmission ratio 32:1
 - > For **endodontic files with shaft Ø 2.35 mm**
- Mat. no. 1.011.6790

The basis for perfect prophylaxis.



SMARTmatic PROPHY S19

- > Transmission ratio 8:1
 - > For attachment of **disposable Prophylaxis heads or steel Prophylaxis heads** with Doriot connection
- Mat. no. 1.011.6740




PROPHYwiz 181 P

- > Prophylactic handpiece with **integrated air motor**
 - > Maximum speed **3,000 min⁻¹**
 - > For attachment of **disposable Prophylaxis heads or steel Prophylaxis heads** with Doriot connection
 - > **Rotates 360°**
- Mat. no. 1.003.2278



SMARTmatic PROPHY S31

- > Transmission ratio 8:1
 - > **Double-threaded (1-64-UNC and 1-72-UNF) for screw-type polisher / brushes.**
 - > **Unscrewable holder** for Snap-On polishing applications
 - > **V-sealing** inside the head for prevention of intrusion of paste
- Mat. no. 1.011.6760



SONICflex 2000 N

- > Airscaler with a frequency of **6,000 Hz**
 - > **Amplitude of 160 µm** for gentle and effective removal of debris through elliptical motions
 - > Selection of more than **50 tips** for a broad range of applications
- Mat. no. 0.571.5020



SMARTmatic PROPHY S33

- > Transmission ratio 8:1
 - > **70° oscillating** Rotation direction
 - > **Holder for Snap-On polishing applications**
 - > **V-sealing** inside the head for prevention of intrusion of paste
- Mat. no. 1.011.6800



SMARTmatic S53

- > Transmission ratio 8:1
 - > **70° oscillating** Rotation direction
 - > For **polishing and root canal instruments with 2,35mm**
- Mat. no. 1.011.6810

12 months warranty.



Look after your investment with perfect handpiece care.



QUATTROcare PLUS

- > **Cost-efficient care for your handpieces**
- > **Optimum quality of care** with the KaVo QUATTROcare plus care spray
- > **Easy efficient handling** and fast programme sequence
- > **Low purchase and consumables costs**

Mat. no. 1.008.3805



KaVo Spray

- > Manual handpiece/ turbine care
- > **Effectively preserves value**

Mat. no. 0.411.9640



You know what you want. We have got what you need.

The KaVo handpiece series –
the perfect handpiece for every dentist.

KaVo master series

Enjoy maximum power with innovative KaVo comfort.

➤ Highest treatment quality

➤ Maximum treatment comfort

➤ Highest patient satisfaction

KaVo expert series

Rely on advanced KaVo performance.

➤ Best durability at unbeatable prices

➤ Proven configuration solutions

➤ Highest efficiency

KaVo smart series

Start safe with affordable KaVo quality.

➤ Made in Germany

➤ Premium materials

➤ Entry level price for the world of KaVo



KaVo. Dental Excellence.